Prior to use, carefully read the instructions for use.

EN

CERASMART™

FORCE ABSORBING HYBRID CERAMIC CAD/CAM BLOCK

For use only by a dental professional in the recommended indications.

INDICATIONS

- Metal free indirect restorations: full crown, inlay, onlay, laminated veneer
- Implant supported crown


CONTRAINDICATIONS

- Avoid use of this product in patients with known allergies to methacrylate polymer.

PREPARATION DESIGN

When designing restorations, the following minimum dimensions should be maintained:

- Minimum wall thickness should be 1.5 mm, 1.0 mm at margins.
- Prepare margins with deep chamfer or rounded shoulder. (Figure 1)
- For inlays and onlays, all internal edges and angles should be rounded. Avoid having margins in direct occlusal contact with the opposing tooth.
- Minimum thickness of the restoration should be 1.5 mm in pit and fissure areas and 1.5 mm in cusp areas.

MILLING

Please refer to respective instructions of the milling unit.

Select the milling program for CERASMART (GC) or composite resin block materials.

After milling, check the restoration for:

- Discolorations linked to the milling process
- Protuberances
- Cracks or chips

If the restoration is defective, discontinue the manufacture of restoration.

FINISHING AND POLISHING

Finish and polish with coarse and fine silicone point. DIAPOLISHER PASTE (GC) can be used for final polishing to obtain high gloss.

NOTE:

When adjusting or polishing the restorative, use a dust collector and wear a dust mask to avoid inhaling dust. Wear protective glasses.

CEMENTATION

With Sandblasting technique

- Sandblasting with 25-50µm alumina 0.15MPa/1.5bar is recommended.
- Blow the restoration with an oil-free air syringe or clean with an ultrasonic cleaner and dry.
- Clean with alcohol to remove oil residue.
- Treat the surface with a silane coupling agent such as CERAMIC PRIMER II (GC).
- Cement with an adhesive resin cement.

NOTE:

Prior to using CERAMIC PRIMER II (GC) or resin cement, refer to the respective instructions for use of the manufacturer.

Without Sandblasting technique

- Treat with hydrofluoric acid (5% for 60 seconds.
- Wash with water spray or an ultrasonic cleaner and dry.
- Treat the surface with a silane coupling agent such as CERAMIC PRIMER II (GC).
- Cement with an adhesive resin cement.

NOTE:

- Prior to using CERAMIC PRIMER II (GC) or resin cement, refer to the respective instructions for use of the manufacturer.
- Hydrofluoric acid is for extraoral use only and must not be applied in the oral cavity.
- Phosphoric acid (35-37%) can be also used for the purpose of cleaning the surface, preferably scrubbing with a microbrush for 10-15 seconds.

ADDITIONAL APPLICATION AND REPAIR

Roughen the bonding area using diamond point or carbide bur. Apply CERAMIC PRIMER II (GC) to improve adhesion and dry the surface. Apply a bonding agent following its manufacturer’s instruction. Proceed with the application of the chosen composite restorative material.

G-aenial Universal Flo (GC) is recommended in order to make the surface of restorative homogeneous as it consists of the same flexible nano ceramic technology as CERASMART. In this case, the intermediate use of a bonding agent is not necessary. Alternatively, OPTIGLAZE Color (GC) is also recommended as surface characterization and glazing technique. Sandblast with 25-50µ alumina (0.15MPa/1.5bar), clean and dry. Apply CERAMIC PRIMER II (GC), dry and apply OPTIGLAZE Color (GC). OPTIGLAZE Color (GC) can only be cured with a light-curing device with a wavelength of 400-430nm. Refer to the respective instructions for use.

SHADE

5 Vita shades* available in two translucencies:

- High Translucency (HT): A1 HT, A2 HT, A3 HT, A3.5 HT, B1 HT
- Low Translucency (LT): A1 LT, A2 LT, A3 LT, A3.5 LT, B1 LT
- One bleach shade (BL)

*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

SIZE

Available in 3 sizes: 12, 14, 14L

PACKAGE

5 blocks per box

STORAGE

Recommended for optimal performance, store at room temperature (4-25°C / 39.2-77.0°F) away from direct sunlight and high humidity. (Shelf life : 5 years from date of manufacture)

CAUTION

- In rare cases this product may cause sensitivity in some people. If such reactions are experienced, discontinue the use of product and consult a physician.

The Safety Data Sheets for this product are available at <http://www.gceurope.com/download> or for North and Latin America <http://www.gcamerica.com/downloads> They can also be obtained from your supplier.

Last revised: 06/2015

Rx Only

0086

GC

MANUFACTURED BY GC DENTAL PRODUCTS CORP. 3737 Torimatsu-cho, Kasugai, Aichi 466-0844, Japan

RESPONSIBLE MANUFACTURER IN CANADA: GC AMERICA INC. 3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED BY GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V. Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL : +32 16 74 10 00

GC AMERICA INC. 3737 West 127th Street, Alsip, IL 60803 U.S.A.

GC ASIA DENTAL PTE. LTD. 11 Tampines Concourse, #03-05, Singapore 528729 TEL : +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD. 1753 Botany Rd, Banksmeadow, NSW 2019, Australia TEL : +61 2 9316 4499

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CONTRAINDICATIONS

1. Avoid use of this product in patients with known allergies to methacrylate polymer.

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- Minimum wall thickness should be 1.5 mm, 1.0 mm at margins. Prepare margins with deep chamfer or rounded shoulder. (Figure 1)
- For inlays and onlays, all internal edges and angles should be rounded. Avoid having margins in direct occlusal contact with the opposing tooth.
- Minimum thickness of the restoration should be 1.5 mm in pit and fissure areas and 1.5 mm in cusps areas.

MILLING

Please refer to respective instructions of the milling unit. Select the milling program for CERASMART (GC) or composite restorative block materials. After milling, check the restoration for:

- Discolorations linked to the milling process
 - Protuberances
 - Cracks or chips
- If the restoration is defective, discontinue the manufacture of restoration.

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NOTE:
When adjusting or polishing the restorative, use a dust collector and wear a dust mask to avoid inhaling dust. Wear protective glasses.

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NOTE:
When using CERAMIC PRIMER II (GC) or resin cement, refer to the respective instructions for use of the manufacturer.

Without Sandblasting technique

- Treat with hydrofluoric acid (5%) for 60 seconds.
 - Wash with water spray or an ultrasonic cleaner and dry.
 - Prict the surface with a silane coupling agent such as CERAMIC PRIMER II (GC).
 - Cement with an adhesive resin cement.
- NOTE:**
- Prior to using CERAMIC PRIMER II (GC) or resin cement, refer to the respective instructions for use of the manufacturer.
 - Hydrofluoric acid should be used only and must not be applied in the oral cavity.
 - Phosphoric acid (35-37%) can be also used for the purpose of cleaning the surface, preferably scrubbing with a microbrush for 10-15 seconds.

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Roughen the bonding area using diamond point or carbide bur. Apply CERAMIC PRIMER II (GC) to improve adhesion and dry the surface. Apply a bonding agent following the manufacturer's instruction. Proceed with the application of the chosen composite restorative material. G-aenial Universal Flo (GC) is recommended in order to make the surface of restorative homogeneous as it consists of the same flexible nano ceramic technology as CERASMART. In this case, the intermediate use of a bonding agent is not necessary. Alternatively, OPTIGLAZE Color (GC) is also recommended as surface characterization and repair. High Translucency (HT): A1 HT, A2 HT, A3 HT, A3.5 HT, B1 HT Low Translucency (LT): A1 LT, A2 LT, A3 LT, A3.5 LT, B1 LT One bleach shade (BL) *Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen,Germany.

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MANUFACTURED BY GC DENTAL PRODUCTS CORP.

DISTRIBUTED BY GC CORPORATION

76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.

2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

PRINTED IN JAPAN

0086

IGI

MANUFACTURED BY GC DENTAL PRODUCTS CORP.

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5 Vita shades* available in two translucencies: - High Translucency (HT): A1 HT, A2 HT, A3 HT, A3.5 HT, B1 HT - Low Translucency (LT): A1 LT, A2 LT, A3 LT, A3.5 LT, B1 LT One bleach shade (BL) *Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen,Germany.